

## NPY antibody

### Product Information

Catalog No.:	FNab05829
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


NPY is implicated in the control of feeding and in secretion of gonadotrophin-release hormone.

### Immunogen information

Immunogen:	neuropeptide Y
Synonyms:	Pro-neuropeptide Y Neuropeptide Y Alternative names: Neuropeptide tyrosine (NPY) C-flanking peptide of NPY (CPON) NPY
Observed MW:	11 kDa
Uniprot ID :	P01303

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab05829 (Neuropeptide Y antibody) at dilution of 1:200