

## NRAS antibody

### Product Information

Catalog No.:	FNab05849
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Ras proteins bind GDP/GTP and possess intrinsic GTPase activity.

### Immunogen information


Immunogen:	neuroblastoma RAS viral(v-ras) oncogene homolog
Synonyms:	GTPase NRas Transforming protein N-Ras NRAS HRAS1
Observed MW:	21 kDa
Uniprot ID :	P01111

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab05849(NRAS Antibody) at dilution of 1:200



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab05849(NRAS Antibody) at dilution of 1:1500