

NRBF2 antibody

Product Information

Catalog No.:	FNab05850
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


Background

Immunogen information


Immunogen:	nuclear receptor binding factor 2
Synonyms:	Nuclear receptor-binding factor 2 (NRBF-2) Comodulator of PPAR and RXR NRBF2 COPR
Observed MW:	38 kDa
Uniprot ID :	Q96F24

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer using FNab05850(NRBF2 Antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05850(NRBF2 Antibody) at dilution of 1:1000