

NRIP3 antibody

Product Information

Catalog No.:	FNab05858
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	nuclear receptor interacting protein 3
Synonyms:	Nuclear receptor-interacting protein 3 Sarcoma antigen NY-SAR-105 NRIP3 C11orf14
Observed MW:	27 kDa
Uniprot ID :	Q9NQ35

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tissue slide using FNab05858(NRIP3 Antibody) at dilution of 1:50



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab05858(NRIP3 Antibody) at dilution of 1:300