

HMGN5 antibody

Product Information

Catalog No.:	FNab05862
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	nucleosomal binding protein 1
Synonyms:	High mobility group nucleosome-binding domain-containing protein 5 Nucleosome-binding protein 1 HMGN5 NSBP1
Observed MW:	55 kDa
Uniprot ID :	P82970

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab05862(NSBP1 Antibody) at dilution of 1:200



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab05862(NSBP1 Antibody) at dilution of 1:600