

ENO2 antibody

Product Information

Catalog No.:	FNab05867
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	0F8
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	enolase 2(gamma, neuronal)
Synonyms:	Gamma-enolase 2-phospho-D-glycerate hydro-lyase Enolase 2 Neural enolase Neuron-specific enolase (NSE) ENO2
Observed MW:	47 kDa
Uniprot ID :	P09104

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:2000-1:20000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab05867(ENO2 Antibody) at dilution of 1:200



U-251 cells were subjected to SDS PAGE followed by western blot with FNab05867(ENO2 Antibody) at dilution of 1:4000