

## NSUN5 antibody

### Product Information

Catalog No.:	FNab05871
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


S-adenosyl-L-methionine-dependent methyltransferase that specifically methylates the C(5) position of cytosine 3782 in 28S rRNA.

### Immunogen information


Immunogen:	NOL1/NOP2/Sun domain family, member 5
Synonyms:	28S rRNA (cytosine-C(5))-methyltransferase NOL1-related protein (NOL1R) NOL1/NOP2/Sun domain family member 5 Williams-Beuren syndrome chromosomal region 20A protein NSUN5 NSUN5A WBSCR20 WBSCR20A
Observed MW:	47 kDa
Uniprot ID :	Q96P11

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta using FNab05871(NSUN5 antibody) at dilution of 1:50



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab05871(NSUN5 antibody) at dilution of 1:500