

## NSUN6 antibody

### Product Information

Catalog No.:	FNab05872
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May have S-adenosyl-L-methionine-dependent methyl-transferase activity.

### Immunogen information

Immunogen:	NOL1/NOP2/Sun domain family, member 6
Synonyms:	tRNA (cytosine(72)-C(5))-methyltransferase NSUN6 NOL1/NOP2/Sun and PUA domain-containing protein 1 NOL1/NOP2/Sun domain family member 6 NSUN6 NOPD1
Observed MW:	52-55 kDa
Uniprot ID :	Q8TEA1

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Neuro-2a cells were subjected to SDS PAGE followed by western blot with FNab05872( NSUN6 Antibody) at dilution of 1:1000