

NT5C3B antibody

Product Information

Catalog No.:	FNab05876
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	5'-nucleotidase, cytosolic III-like
Synonyms:	7-methylguanosine phosphate-specific 5'-nucleotidase (7-methylguanosine nucleotidase) Cytosolic 5'-nucleotidase 3B Cytosolic 5'-nucleotidase III-like protein (cN-III-like protein) N(7)-methylguanylate 5'-phosphatase NT5C3B NT5C3L
Observed MW:	30-34 kDa
Uniprot ID :	Q969T7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



human testis tissue were subjected to SDS PAGE followed by western blot with FNab05876(NT5C3L antibody) at dilution of 1:800