

NTM antibody

Product Information

Catalog No.:	FNab05885
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Neural cell adhesion molecule.

Immunogen information

Immunogen:	neurotrimin
Synonyms:	Neurotrimin (hNT) IgLON family member 2 NTM IGLON2 NT
Observed MW:	58 kDa, 65 kDa
Uniprot ID :	Q9P121

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab05885(NTM Antibody) at dilution of 1:600