

GNL3 antibody

Product Information

Catalog No.:	FNab05894
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	guanine nucleotide binding protein-like 3(nucleolar)
Synonyms:	Guanine nucleotide-binding protein-like 3 E2-induced gene 3 protein Novel nucleolar protein 47 (NNP47) Nucleolar GTP-binding protein 3 Nucleostemin GNL3 E2IG3 NS
Observed MW:	60-62 kDa
Uniprot ID :	Q9BVP2

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab05894(Nucleostemin antibody) at dilution of 1:2000