

NUDT19 antibody

Product Information

Catalog No.:	FNab05902
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Coenzyme A diphosphatase that mediates the hydrolysis of a wide range of CoA esters, including choloyl-CoA and branched-chain fatty-acyl-CoA esters. At low substrate concentrations medium and long-chain fatty-acyl-CoA esters are the primary substrates(By similarity).

Immunogen information


Immunogen:	nudix(nucleoside diphosphate linked moiety X)-type motif 19
Synonyms:	Acyl-coenzyme A diphosphatase NUDT19 Nucleoside diphosphate-linked moiety X motif 19 (Nudix motif 19) NUDT19
Observed MW:	42-45 kDa
Uniprot ID :	A8MXV4

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab05902(NUDT19 antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab05902(NUDT19 antibody) at dilution of 1:1500