

## NUDT8 antibody

### Product Information

Catalog No.:	FNab05907
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Probably mediates the hydrolysis of some nucleoside diphosphate derivatives.

### Immunogen information


Immunogen:	nudix(nucleoside diphosphate linked moiety X)-type motif 8
Synonyms:	Mitochondrial coenzyme A diphosphatase NUDT8 Nucleoside diphosphate-linked moiety X motif 8 (Nudix motif 8) NUDT8
Observed MW:	25 kDa
Uniprot ID :	Q8WV74

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab05907(NUDT8 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab05907(NUDT8 antibody) at dilution of 1:1500