

## NUFIP2 antibody

### Product Information

Catalog No.:	FNab05911
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Binds RNA.

### Immunogen information


Immunogen:	nuclear fragile X mental retardation protein interacting protein 2
Synonyms:	FMR1-interacting protein NUFIP2 82 kDa FMRP-interacting protein (82-FIP) Cell proliferation-inducing gene 1 protein FMRP-interacting protein 2 Nuclear FMR1-interacting protein 2 NUFIP2 KIAA1321
Observed MW:	76-82 kDa
Uniprot ID :	Q7Z417

### Application

Reactivity:	Human
Tested Application:	ELISA, IP, WB
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000
Image:	



IP Result of anti-NUFIP2 (IP:FNab05911, 4ug;  
Detection:FNab05911 1:500) with HeLa cells  
lysate 2000ug.



HeLa cells were subjected to SDS PAGE followed  
by western blot with FNab05911(NUFIP2  
Antibody) at dilution of 1:600