

NUMA1 antibody

Product Information

Catalog No.:	FNab05912
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	nuclear mitotic apparatus protein 1
Synonyms:	Nuclear mitotic apparatus protein 1 Nuclear matrix protein-22 (NMP-22) Nuclear mitotic apparatus protein (NuMA protein) SP-H antigen NUMA1 NMP22 NUMA
Observed MW:	238 kDa
Uniprot ID :	Q14980

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab05912(NUMA antibody) at dilution of 1:100



Hela cells were subjected to SDS PAGE followed by western blot with FNab05912(NUMA antibody) at dilution of 1:1000