

NUP160 antibody

Product Information

Catalog No.:	FNab05922
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in poly(A)+ RNA transport.

Immunogen information

Immunogen:	nucleoporin 160kDa
Synonyms:	Nuclear pore complex protein Nup160 160 kDa nucleoporin Nucleoporin Nup160 NUP160 KIAA0197 NUP120
Observed MW:	149 kDa
Uniprot ID :	Q12769

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05922(NUP160 antibody) at dilution of 1:500