

NUP205 antibody

Product Information

Catalog No.:	FNab05923
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a role in the nuclear pore complex(NPC) assembly and/or maintenance. May anchor NUP62 and other nucleoporins, but not NUP153 and TPR, to the NPC.

Immunogen information


Immunogen:	nucleoporin 205kDa
Synonyms:	Nuclear pore complex protein Nup205 205 kDa nucleoporin Nucleoporin Nup205 NUP205 C7orf14 KIAA0225
Observed MW:	205 kDa
Uniprot ID :	Q92621

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney slide using FNab05923(NUP205 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05923(NUP205 Antibody) at dilution of 1:2000