

NUP54 antibody

Product Information

Catalog No.:	FNab05926
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Component of the nuclear pore complex, a complex required for the trafficking across the nuclear membrane.

Immunogen information

Immunogen:	nucleoporin 54kDa
Synonyms:	Nucleoporin p54 54 kDa nucleoporin NUP54
Observed MW:	55 kDa
Uniprot ID :	Q7Z3B4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05926(NUP54 antibody) at dilution of 1:2000