

## NUP58 antibody

### Product Information

|               |  |
|---------------|--|
| Catalog No.:  | FNab05933  |
| Size:         | 100µg  |
| Form:         | liquid   |
| Purification: | Immunogen affinity purified  |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Rabbit   |
| Clonality:    | polyclonal   |
| Clone ID:     | None   |
| IsoType:      | IgG  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Background


None

### Immunogen information

|              |  |
|--------------|--|
| Immunogen:   | nucleoporin like 1   |
| Synonyms:    | Nucleoporin p58/p45 58 kDa nucleoporin Nucleoporin-like protein 1 NUP58 KIAA0410 NUPL1 |
| Observed MW: | 65 kDa   |
| Uniprot ID : | Q9BVL2   |

### Application

|                       |                   |
|-----------------------|-------------------|
| Reactivity:           | Human, Mouse, Rat |
| Tested Application:   | ELISA, WB         |
| Recommended dilution: | WB: 1:500-1:2000  |
| Image:                |                   |



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab05933(NUPL1 antibody) at dilution of 1:500