

NR4A1 antibody

Product Information

Catalog No.:	FNab05935
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	nuclear receptor subfamily 4, group A, member 1
Synonyms:	Nuclear receptor subfamily 4immunitygroup A member 1 Early response protein NAK1 Nuclear hormone receptor NUR/77 (Nur77) Orphan nuclear receptor HMR Orphan nuclear receptor TR3 ST-59 Testicular receptor 3 NR4A1 GFRP1 HMR NAK1
Observed MW:	70-80 kDa
Uniprot ID :	P22736

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab05935(NUR77 antibody) at dilution of 1:600