

## NVL antibody

### Product Information

Catalog No.:	FNab05939
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Participates in the assembly of the telomerase holoenzyme and effecting of telomerase activity via its interaction with TERT(PubMed:22226966). May play a role in 60S ribosomal subunit biogenesis(PubMed:15469983, PubMed:16782053).

### Immunogen information

Immunogen:	nuclear VCP-like
Synonyms:	Nuclear valosin-containing protein-like (NVLp, Nuclear VCP-like protein) NVL NVL2
Observed MW:	95 kDa
Uniprot ID :	O15381

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse thymus tissue were subjected to SDS PAGE followed by western blot with FNab05939(NVL antibody) at dilution of 1:3000