

## NWD1 antibody

### Product Information

Catalog No.:	FNab05940
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in the control of androgen receptor(AR) protein steady-state levels.

### Immunogen information

Immunogen:	NACHT and WD repeat domain containing 1
Synonyms:	NACHT domain- and WD repeat-containing protein 1 NWD1
Observed MW:	77 kDa
Uniprot ID :	Q149M9

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:1000; IF: 1:20-1:200
Image:	



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab05940( NWD1 Antibody) at dilution of 1:600