

## NYAP2 antibody

### Product Information

Catalog No.:	FNab05947
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Activates PI3K and concomitantly recruits the WAVE1 complex to the close vicinity of PI3K and regulates neuronal morphogenesis.

### Immunogen information


Immunogen:	KIAA1486 protein
Synonyms:	Neuronal tyrosine-phosphorylated phosphoinositide-3-kinase adapter 2 NYAP2 KIAA1486
Observed MW:	80 kDa
Uniprot ID :	Q9P242

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab05947(KIAA1486 Antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab05947(KIAA1486 Antibody) at dilution of 1:1000