

POU2F2-Specific antibody

Product Information

Catalog No.:	FNab05962
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	POU class 2 homeobox 2
Synonyms:	POU domain, class 2, transcription factor 2 Lymphoid-restricted immunoglobulin octamer-binding protein NF-A2 Octamer-binding protein 2 (Oct-2) Octamer-binding transcription factor 2 (OTF-2) POU2F2 OCT2 OTF2
Observed MW:	55 kDa
Uniprot ID :	P09086

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human spleen using FNab05962(OCT2-Specific antibody) at dilution of 1:100



mouse thymus tissue were subjected to SDS PAGE followed by western blot with FNab05962(OCT2-Specific antibody) at dilution of 1:300