

POU5F1 antibody

Product Information

Catalog No.:	FNab05964
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	0H7
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	POU class 5 homeobox 1
Synonyms:	POU domain, class 5, transcription factor 1 Octamer-binding protein 3 (Oct-3) Octamer-binding protein 4 (Oct-4) Octamer-binding transcription factor 3 (OTF-3) POU5F1 OCT3 OCT4 OTF3
Observed MW:	50 kDa
Uniprot ID :	Q01860

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:4000-1:20000; IHC: 1:200-1:1000
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab05964(OCT4 Antibody) at dilution of 1:500



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab05964(OCT4 Antibody) at dilution of 1:10000