

ODC1 antibody

Product Information

Catalog No.:	FNab05968
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Key enzyme of polyamine biosynthesis that converts ornithine into putrescine, which is the precursor for the polyamines, spermidine and spermine.

Immunogen information


Immunogen:	ornithine decarboxylase 1
Synonyms:	Ornithine decarboxylase (ODC) ODC1
Observed MW:	53 kDa
Uniprot ID :	P11926

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate hyperplasia using FNab05968(ODC1 antibody) at dilution of 1:50



mouse thymus tissue were subjected to SDS PAGE followed by western blot with FNab05968(ODC1 Antibody) at dilution of 1:300