

ODF2L antibody

Product Information

Catalog No.:	FNab05972
Size:	100μg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	outer dense fiber of sperm tails 2-like
Synonyms:	Protein BCAP Basal body centriole-associated protein Outer dense fiber protein 2-like (ODF2-like protein) ODF2L BCAP KIAA1229
Observed MW:	60-74 kDa
Uniprot ID :	Q9ULJ1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab05972(ODF2L antibody) at dilution of 1:600