

ZNF830 antibody

Product Information


Catalog No.:	FNab05992
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	zinc finger protein 830
Synonyms:	Zinc finger protein 830 Coiled-coil domain-containing protein 16 ZNF830 CCDC16
Observed MW:	42 kDa
Uniprot ID :	Q96NB3

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05992(ZNF830) antibody at dilution of 1:1000