

OSBPL9 antibody

Product Information

Catalog No.:	FNab06022
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	oxysterol binding protein-like 9
Synonyms:	Oxysterol-binding protein-related protein 9 (ORP-9, OSBP-related protein 9) OSBPL9 ORP9 OSBP4
Observed MW:	90 kDa
Uniprot ID :	Q96SU4

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab06022(OSBPL9 antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab06022(OSBPL9 antibody) at dilution of 1:300