

ATP5PO antibody

Product Information

Catalog No.:	FNab06023
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ATP synthase, H ⁺ transporting, mitochondrial F1 complex, O subunit
Synonyms:	ATP synthase subunit O, mitochondrial ATP synthase peripheral stalk subunit OSCP Oligomycin sensitivity conferral protein (OSCP) ATP5PO ATP5O ATPO
Observed MW:	23-25 kDa
Uniprot ID :	P48047

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab06023(ATP5O antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab06023(ATP5O antibody) at dilution of 1:1000