

## BGLAP antibody

### Product Information

Catalog No.:	FNab06030
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	bone gamma carboxyglutamate protein 1
Synonyms:	Osteocalcin Bone Gla protein (BGP) Gamma-carboxyglutamic acid-containing protein BGLAP
Observed MW:	None
Uniprot ID :	P02818

### Application

Reactivity:	Mouse
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded mouse testis tissue using FNab06030(Osteocalcin antibody) at dilution of 1:200