

OSTM1 antibody

Product Information

Catalog No.:	FNab06036
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Required for osteoclast and melanocyte maturation and function.

Immunogen information


Immunogen:	osteopetrosis associated transmembrane protein 1
Synonyms:	Osteopetrosis-associated transmembrane protein 1 Chloride channel 7 beta subunit OSTM1 GL
Observed MW:	35-45 kDa
Uniprot ID :	Q86WC4

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-OSTM1 (IP:FNab06036, 4ug;
Detection:FNab06036 1:500) with mouse brain
tissue lysate 4000ug.



mouse brain tissue were subjected to SDS PAGE
followed by western blot with FNab06036(OSTM1
antibody) at dilution of 1:1000