

OTUD4 antibody

Product Information

Catalog No.:	FNab06040
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Deubiquitinating enzyme that specifically hydrolyzes 'Lys-48'-linked polyubiquitin.

Immunogen information


Immunogen:	OTU domain containing 4
Synonyms:	OTU domain-containing protein 4 HIV-1-induced protein HIN-1 OTUD4 HIN-1 KIAA1046
Observed MW:	140-150 kDa
Uniprot ID :	Q01804

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab06040(OTUD4 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06040(OTUD4 Antibody) at dilution of 1:1000