

OTUD6A antibody

Product Information

Catalog No.:	FNab06041
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Deubiquitinating enzyme that hydrolyzes 'Lys-27'-, 'Lys-29'-and 'Lys-33'-linked polyubiquitin chains. Also able to hydrolyze 'Lys-11'-linked ubiquitin chains.

Immunogen information

Immunogen:	OTU domain containing 6A
Synonyms:	OTU domain-containing protein 6A DUBA-2 OTUD6A DUBA2
Observed MW:	33 kDa
Uniprot ID :	Q7L8S5

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab06041(OTUD6A Antibody) at dilution of 1:1000