

CTNND1 antibody

Product Information

Catalog No.:	FNab06060
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3F8
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	catenin, delta 1
Synonyms:	Catenin delta-1 Cadherin-associated Src substrate (CAS) p120 catenin (p120(ctn)) p120(cas) CTNND1 KIAA0384
Observed MW:	120 kDa
Uniprot ID :	O60716

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:100-1:400
Image:	



Immunohistochemistry of paraffin-embedded human breast ductal carcinoma tissue slide using FNab06060(p120 Catenin antibody) at dilution of 1:200



fetal human brain tissue were subjected to SDS PAGE followed by western blot with FNab06060(p120 Catenin antibody at dilution of 1:1000