

DCTN1 antibody

Product Information

Catalog No.:	FNab06062
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	dynactin 1(p150, glued homolog, Drosophila)
Synonyms:	Dynactin subunit 1 150 kDa dynein-associated polypeptide DAP-150 (DP-150) p135 p150-glued DCTN1
Observed MW:	150 kDa
Uniprot ID :	Q14203

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab06062(DCTN1 antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab06062(DCTN1 antibody) at dilution of 1:500