

RPRD1A antibody

Product Information

Catalog No.:	FNab06063
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	regulation of nuclear pre-mRNA domain containing 1A
Synonyms:	Regulation of nuclear pre-mRNA domain-containing protein 1A Cyclin-dependent kinase inhibitor 2B-related protein p15INK4B-related protein RPRD1A P15RS
Observed MW:	36 kDa
Uniprot ID :	Q96P16

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab06063(P15RS antibody) at dilution of 1:1000