

AMACR antibody

Product Information

Catalog No.:	FNab06082
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2C4
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	alpha-methylacyl-CoA racemase
Synonyms:	Alpha-methylacyl-CoA racemase 2-methylacyl-CoA racemase AMACR
Observed MW:	42 kDa
Uniprot ID :	Q9UHK6

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab06082(p504S,AMACR antibody) at dilution of 1:1000