

PACS2-Specific antibody

Product Information

Catalog No.:	FNab06099
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	phosphofurin acidic cluster sorting protein 2
Synonyms:	Phosphofurin acidic cluster sorting protein 2 (PACS-2) PACS1-like protein PACS2 KIAA0602 PACS1L
Observed MW:	116 kDa
Uniprot ID :	Q86VP3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human liver tissue slide using FNab06099(PACS2-Specific antibody) at dilution of 1:200



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab06099(PACS2-Specific antibody) at dilution of 1:600