

## PAFAH1B3 antibody

### Product Information

Catalog No.:	FNab06109
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

Inactivates paf by removing the acetyl group at the sn-2 position. This is a catalytic subunit. Plays an important role during the development of brain.

### Immunogen information

Immunogen:	platelet-activating factor acetylhydrolase, isoform Ib, gamma subunit 29kDa
Synonyms:	Platelet-activating factor acetylhydrolase IB subunit alpha1 PAF acetylhydrolase 29 kDa subunit (PAF-AH 29 kDa subunit) PAF-AH subunit gamma (PAFAH subunit gamma) PAFAH1B3 PAFAHG
Observed MW:	None
Uniprot ID :	Q15102

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA