

PALLD-Specific antibody

Product Information

Catalog No.:	FNab06127
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	palladin, cytoskeletal associated protein
Synonyms:	Palladin SIH002 Sarcoma antigen NY-SAR-77 PALLD KIAA0992
Observed MW:	120 kDa
Uniprot ID :	Q8WX93

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:1000; IP: 1:200-1:1000
Image:	



IP Result of anti-Palladin-Specific (IP:FNab06127, 4ug; Detection:FNab06127 1:200) with HeLa cells lysate 1500ug.



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06127(Palladin-Specific antibody) at dilution of 1:200