

## PAMR1 antibody

### Product Information

Catalog No.:	FNab06130
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in regeneration of skeletal muscle.

### Immunogen information


Immunogen:	regeneration associated muscle protease
Synonyms:	Inactive serine protease PAMR1 Peptidase domain-containing protein associated with muscle regeneration 1 Regeneration-associated muscle protease homolog PAMR1 RAMP
Observed MW:	80 kDa
Uniprot ID :	Q6UXH9

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle using FNab06130(PAMR1 antibody) at dilution of 1:50



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab06130(PAMR1 antibody) at dilution of 1:500