

## PCSK5 antibody

### Product Information

Catalog No.:	FNab06186
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	proprotein convertase subtilisin/kexin type 5
Synonyms:	Proprotein convertase subtilisin/kexin type 5 Proprotein convertase 5 (PC5) Proprotein convertase 6 (PC6, hPC6) Subtilisin/kexin-like protease PC5 PCSK5 PC5 PC6
Observed MW:	202 kDa
Uniprot ID :	Q92824

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab06186(PC5 antibody) at dilution of 1:200



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab06186(PC5 Antibody) at dilution of 1:300