

KAT2B antibody

Product Information


Catalog No.:	FNab06188
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	K(lysine) acetyltransferase 2B
Synonyms:	Histone acetyltransferase KAT2B Histone acetyltransferase PCAF (Histone acetylase PCAF) Lysine acetyltransferase 2B P300/CBP-associated factor (P/CAF) Spermidine acetyltransferase KAT2B KAT2B PCAF
Observed MW:	95-110 kDa
Uniprot ID :	Q92831

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab06188 at dilution of 1:500