

PCDH11X antibody

Product Information

Catalog No.:	FNab06195
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Potential calcium-dependent cell-adhesion protein.

Immunogen information


Immunogen:	protocadherin 11 X-linked
Synonyms:	Protocadherin-11 X-linked (Protocadherin-11) Protocadherin on the X chromosome (PCDH-X) Protocadherin-S PCDH11X KIAA1326 PCDH11 PCDHX
Observed MW:	135 kDa
Uniprot ID :	Q9BZA7

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab06195(PCDH11X antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab06195(PCDH11X antibody) at dilution of 1:400