

## PCDHA2 antibody

### Product Information

Catalog No.:	FNab06198
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Potential calcium-dependent cell-adhesion protein. May be involved in the establishment and maintenance of specific neuronal connections in the brain.

### Immunogen information


Immunogen:	protocadherin alpha 2
Synonyms:	Protocadherin alpha-2 (PCDH-alpha-2) PCDHA2
Observed MW:	80 kDa
Uniprot ID :	Q9Y5H9

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-PCDHA2 (IP:FNab06198, 3ug;  
Detection:FNab06198 1:300) with HeLa cells  
lysate 4650ug.



human brain tissue were subjected to SDS PAGE  
followed by western blot with  
FNab06198(PCDHA2 antibody) at dilution of  
1:800