

PCDHA3 antibody

Product Information

Catalog No.:	FNab06199
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Potential calcium-dependent cell-adhesion protein. May be involved in the establishment and maintenance of specific neuronal connections in the brain.

Immunogen information


Immunogen:	protocadherin alpha 3
Synonyms:	Protocadherin alpha-3 (PCDH-alpha-3) PCDHA3
Observed MW:	102 kDa
Uniprot ID :	Q9Y5H8

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab06199(PCDHA3 antibody) at dilution of 1:50



A549 cells were subjected to SDS PAGE followed by western blot with FNab06199(PCDHA3 antibody) at dilution of 1:300