

## PCDHB5 antibody

### Product Information

Catalog No.:	FNab06204
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Potential calcium-dependent cell-adhesion protein. May be involved in the establishment and maintenance of specific neuronal connections in the brain.

### Immunogen information

Immunogen:	protocadherin beta 5
Synonyms:	Protocadherin beta-5 (PCDH-beta-5) PCDHB5
Observed MW:	70-75 kDa
Uniprot ID :	Q9Y5E4

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



K562 cells were subjected to SDS PAGE followed by western blot with FNab06204(PCDHB5 antibody) at dilution of 1:800