

PCDHGB3 antibody

Product Information

Catalog No.:	FNab06205
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Potential calcium-dependent cell-adhesion protein. May be involved in the establishment and maintenance of specific neuronal connections in the brain.

Immunogen information

Immunogen:	protocadherin gamma subfamily B, 3
Synonyms:	Protocadherin gamma-B3 (PCDH-gamma-B3) PCDHGB3
Observed MW:	89 kDa
Uniprot ID :	Q9Y5G1

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:1000; IF: 1:10-1:100
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab06205(PCDHGB3 antibody) at dilution of 1:400