

PCF11 antibody

Product Information

Catalog No.:	FNab06207
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Component of pre-mRNA cleavage complex II.

Immunogen information

Immunogen:	PCF11, cleavage and polyadenylation factor subunit, homolog
Synonyms:	Pre-mRNA cleavage complex 2 protein Pcf11 Pre-mRNA cleavage complex II protein Pcf11 PCF11 KIAA0824
Observed MW:	173 kDa
Uniprot ID :	O94913

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



K562 cells were subjected to SDS PAGE followed by western blot with FNab06207(PCF11 antibody) at dilution of 1:800